## DAFTAR PUSTAKA

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### Intended use:

Enzymatic in vitro test for the quantitative determination of uric acid in human serum and plasma.

### Summary:


Uric acid is the final product of purine metabolism in the human organism. Uric acid measurements are used in the diagnosis and treatment of numerous renal and metabolic disorders, including renal failure, gout, leukaemia, psoriasis, starvation or other wasting conditions, and of patients receiving cytotoxic drugs. The oxidation of uric acid provides for two approaches to the quantitative determination of this purine metabolite. One approach is the reduction of phosphotungstic acid in an alkaline solution to tungsten blue, which is measured photometrically. The method is, however, subject to interferences from drugs and reducing substances other than uric acid.

A second approach, described by Praetorius and Poulsen, utilizes the enzyme uricase to oxidize uric acid; this method eliminates the interferences intrinsic to chemical oxidation. Uricase can be employed in methods that involve the UV measurement of the consumption of uric acid or in combination with other enzyme to provide a colorimetric method.

The assay described here is a slight modification of the colorimetric method. The modifications were described by Siedel. This reaction, the peroxide reacts in the presence of peroxidase, DHBSA and 4-Aminonaphtyrene to form a quinoneimine dye. The intensity of the red color is proportional to the uric acid concentration and is determined photometrically.

### Test principle:

Colorimetric endpoint assay (Uricase-PAP-method)

The increase in absorbance is measured.  
Permeation absorbance & kurve.

### Reagent concentration:

R1:	
Phosphate buffer pH 7.4	50 mmol/l
DHBSA*	4 mmol/l
Uricase	60 U/l
POD	660 U/l
4-Aminonaphtyrene	1 mmol/l
Preservative	
*3,5-Dichloro-2-hydroxy-benzenesulfonic acid	

### R4:

Uric acid	6 mg/dl (356.9 µmol/l)

### Preparation and stability:

R1: ready for use

R4: ready for use

The reagents are stable up to the expiry date in the label when stored light protected at +2°C to +8°C.

After opening: to expire date at +2°C to +8°C, protect from light  
21 days at +20°C to +25°C, protect from light

Coloration of the reagent (reagent blank at 546 nm, 1 cm > 0.2) indicates a contamination, damage during transport or due storage at higher temperatures.

### Specimen:

Serum/plasma

Collect serum using standard sampling tubes.

Heparin, or EDTA-plasma

Stability: 5 days at +2°C to +8°C

6 months at -20°C

### Urine

Collect urine without using preservatives.

Stability: Assay urinary uric acid as soon as possible. Do not refrigerate.

Centrifuge samples containing precipitate before performing the assay.

### Notes:

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

### Limitations - interference:

Criterion: Recovery within ±10% of initial values

Icterus: No significant interference up to an index I of 12 (approximate conjugated and unconjugated bilirubin concentration 12 mg/dl)

Hemolysis: No significant interference up to an index H of 50 (approximate hemoglobin concentration 50 mg/dl).

Lipemia (Intralipid). Triglycerides may interfere. There is poor correlation between turbidity and triglycerides concentration.

Elevated levels of ascorbic acid produces false low values.

Uricase reacts specifically with uric acid. Other purine derivatives can inhibit the uric acid reaction.

### Testing procedure:

Applications for automated systems are available on request

#### Materials provided

- Working solutions as described above

#### Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

#### Manual procedure:

Wavelength:	Hg[546 nm] (490 - 550)
Temperature:	+25 / +30 / +37°C
Cuvette:	1 cm light path
Zero adjustment:	against reagent blank one reagent blank per series only

	Blank	Calibrator	Sample
Calibrator/R4	—	20 µl	—
Sample	—	—	20 µl
R1	1000 µl	1000 µl	1000 µl

Mix, incubate for 5 min. at +37°C or 10 min at +20°C or +25°C.  
Read absorbance of the sample against reagent blank within 30 min.

#### Calculation:

By calibrator:

$$\frac{\Delta A_{\text{sample}}}{\Delta A_{\text{calibrator}}} \times \text{calibrator conc.} = \text{Uric acid conc.}$$

Unit conversion:  
mg/dl × 59.5 = µmol/l  
mg/dl × 0.059 = mmol/l

#### Measuring /reportable range:

0.2 – 20.0 mg/dl (11.9 – 1190 µmol/l)

Determine samples with uric acid concentrations > 20 mg/dl via the rerun function. On instruments without rerun function, manually dilute the samples with 0.9% NaCl or distilled/deionized water (e.g. 1 + 1). Multiply the result by the appropriate dilution factor (e.g. factor 2).

#### Expected values:

Serum/plasma:

Male: 3.4 – 7.0 mg/dl (202.3 – 416.5 µmol/l)

Female: 2.4 – 5.7 mg/dl (142.8 – 339.2 µmol/l)

#### Urine

(Reference range according to Krieg and Colombo)

Morning urine 37 – 92 mg/dl (2200 – 5475 µmol/l)

24 hour urine 200 – 1000 mg/24h (1200 – 5900 µmol/l/24h)

corresponding to 13 – 67 mg/dl (773 – 3986 µmol/l\*)

\*Calculated from a urine volume of 1.5 l/24h

Urine (Reference range according to Tietz)

Average diet: 250 – 750 mg/24h

Low purine diet: Male: < 480 mg/24h

Female: < 400 mg/24h

High purine diet: < 1000 mg/24h

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the uric acid results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

#### Analytical sensitivity (lower detection limit):

Detection limit: 0.2 mg/dl (11.9 µmol/l)

The lower detection limit represents the lowest measurable uric acid concentration that can be distinguished from zero.

# Fluitest® UA

URIC ACID

#### Imprecision:

##### Serum

Reproducibility was determined using human samples and controls ( $n = 20$ ). The following results were obtained:

Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Sample1	4.23	0.03	0.74
Sample2	5.96	0.04	0.65
Sample3	11.49	0.06	0.53

Reproducibility was determined using human samples and controls ( $n = 20$ ). The following results were obtained:

Sample	Between day		
	Mean mg/dl	SD mg/dl	CV %
Sample1	4.70	0.04	0.92
Sample2	6.44	0.08	1.27
Sample3	11.03	0.13	1.20

#### Method comparison:

A comparison of the Analyticon Fluitest® UA (y) with a commercial obtainable assay (x) gave with 44 samples the following result:

$$y = 0.897 x - 1.149; \quad r = 0.9885$$

#### Quality Control:

##### Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
2C Contropath® Plus	2C x 5 ml	#1220
	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

#### Calibration:

S1:	0.9% NaCl		
S2:	Bio Cal®	20 x 3 ml	# 1420
	Bio Cal® E	10 x 3 ml	# 1430

R4: calibrator provided in kit

#### Calibration frequency:

Two-point calibration is recommended

- after lot change
- as required following quality control procedures

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Analyticon® Biotechnologies AG
(AB-GB-UA-01/2)