

DAFTAR PUSTAKA

- Badan Pengawas Obat dan Makanan. 2022. *Pedoman Mitigasi Risiko Cemaran Etilen Glikol (EG) dan Dietilen Glikol (DEG) pada Pangan Olahan*. Badan Pengawas Obat dan Makanan Republik Indonesia. Jakarta. Hal. 1, 20 – 22.
- Baffi, P., Elneser, S., Baffi, M., and De Melin, M. 2000. Quantitative Determination of Diethylene Glycol Contamination in Pharmaceutical Products. *Journal of AOAC International*, 83(4): 793 – 801.
- Belouafa, S., Habti, F., Benhar, S., Belafikh, B., Tayane, S., Hamdouch, S., et al. 2017. Statistical Tools and Approaches to Validate Analytical Methods: Methodology and Practical Examples. *International Journal of Metrology and Quality Engineering*. 8. (9): 1 – 10. doi: 10.1051/ijmqe/2016030
- Bhagat, V., Bhosale, S.T., Kengar, M.D., Jadhav, V.S., and Mohite, M.M. 2023. Ethylene Glycol and Diethylene Glycol Responsible for Adulteration of Cough Syrup. *International Journal of Current Science (IJCSUB)*. 13. (1): 754 – 758.
- Caldeira, L.R., Madureira, F.D., Maia, T.D.F., Muller, V.M., and Fernandes, C. 2021. Simultaneous Quantification of Ethylene Glycol and Diethylene Glycol in Beer by Gas Chromatography Coupled to Mass Spectrometry. *Food Chemistry*. 346. doi: 10.1016/j.foodchem.2020.128871
- Carney, E.W., Freshour, N.L., Dittenber, D.A., and Dryzga, D. 1999. Ethylene Glycol Developmental Toxicity: Unraveling the Roles of Glycolic Acid and Metabolic Acidosis. *Toxicological Sciences*. 50: 117 – 126.
- Chen, Y-H., Liu, H-P., Chen, H-Y., Tsai, F-J., Chang, C-H., and Lee, Y-J. 2011. Ethylene Glycol Induces Calcium Oxalate Crsytal Deposition in Malpighian Tubules: A Drosophila Model for Nephrolithiasis/Urolithiasis. *Kidney International*, 80(4): 369 – 377.
- Dachriyanus. 2004. *Analisis Struktur Senyawa secara Spektroskopi*. Lembaga Pengembangan Teknologi Informasi dan Komunikasi (LPTIK) Universitas Andalas. Padang.
- David, V., Moldoveanu, S.C., and Galaon, T. 2020. Derivatization Procedures and Their Analytical Performances for HPLC Determination in Bioanalysis. doi: 0.1002/bmc.5008
- Decroix, G.A., Gobert, J.G., and DeDeurwaerder, R. 1968. Gas Chromatographic Method for The Determination of Glycerol in Incubates of Adipose Tissue. *Anal. Biochem*. 25: 523 – 531.
- Depkes RI. 2020. *Farmakope Indonesia Edisi VI*. Kementerian Kesehatan Republik Indonesia. Jakarta.
- Dziadosz, M. 2018. Direct Analysis of Ethylene Glycol in Human Serum on the Basis of Analyte Adduct Formation and Liquid Chromatography-Tandem Mass Spectrometry. *Journal of Chromatography B*. 1072: 100 – 104. doi: 10.1016/j.jchromb.2017.11.004

- Ershadi, S., and Shayanfar, A. 2018. Are LOD and LOQ Reliable Parameters for Sensitivity Evaluation of Spectroscopic Methods?. *Journal of AOAC International*. 101. (4): 1213. doi: 10.5740/jaoacint.17-0363
- Fanali, S., Chankvetadze, B., Haddad, P.R., Poole, C.F., and Riekkola, M.L. 2023. *Handbooks in Separation Science Liquid Chromatography: Fundamentals and Instrumentation*. Vol 1, 3rd ed. Elsevier, Belanda.
- Ferenczi-Fodor, K., Vegh, Z., and Nagy-Turak, A. 2001. Validation and Quality Assurance of Planar Chromatographic Procedures in Pharmaceutical Analysis. *Journal of AOAC International*. 84. (4): 1265 – 1276.
- Fleitas-Salazar, N., Silva-Campa, E., Pedroso-Santana, S., Tanori, J., Pedroza-Montero, M.R., and Riera, R. 2017. Effect of Temperature on The Synthesis of Silver Nanoparticles with Polyethylene Glycol: New Insights into The Reduction Mechanism. *J Nanopart Res*, 19(113): 1 – 12. doi: 10.1007/s11051-017-3780-3
- Fowles, J., Banton, M., Klapacz, J., and Shen, H. 2017. A Toxicological Review of The Ethylene Glycol Series: Commonalities and Differences in Toxicity and Modes of Action. *Toxicology Letters*, 278: 66–83. doi: 10.1016/j.toxlet.2017.06.009
- Fu, Y., Hao, Z., Parker, B., and Knapp, M. 2011. A Greener Way to Screen Toothpaste for Diethylene Glycol. *American Journal of Analytical Chemistry*. 2: 938 – 943. doi: 10.1007/s00216-008-2082-0
- Gad, S.C. 2008. *Pharmaceutical Manufacturing Handbook Production and Processes*. John Wiley & Sons, Inc., Hoboken, New Jersey.
- Galli, V., and Barbas, C. 2004. High-Performance Liquid Chromatographic Analysis of Dextromethorphan, Guaifenesin and Benzoate in A Cough Syrup for Stability Testing. *Journal of Chromatography A*, 1048: 207 – 211. doi: 10.1016/j.chroma.2004.07.050
- Ghai, A., Behnke, B., Schulte, K., Jones, T., Goker, D., and Scherbak, D. 2022. An Acute Poisoning of Ethylene Glycol. *EC Clinical and Medical Case Reports*. 5(6): 65 – 68.
- Helawa, I.L. 1994. Stability – Indicating HPLC Assay for Paracetamol, Guaifenesin, Sodium Benzoate and Oxememazine in Cough Syrup. *Analytical Letters*, 27(1): 71 – 93. doi: 10.1080/00032719408006347
- Hernández, F., Ibáñez, M., and Sancho, J.V. 2008. Fast Determination of Toxic Diethylene Glycol in Toothpaste by Ultra-Performance Liquid Chromatography–Time of Flight Mass Spectrometry. *Anal Bioanal Chem*. 391: 1021 – 1027. doi: 10.1007/s00216-008-2082-0
- Holloway, G., Maheswaran, R., Leeks, A., Bradby, S., and Wahab, S. 2010. Screening Method for Ethylene Glycol and Diethylene Glycol in Glycerin-Containing Products. *Journal of Pharmaceutical and Biomedical Analysis*. 51: 507 – 511.

- ICH, (International Conference on Harmonization). 2023. Validation of Analytical Procedures: ICH Guidelines Q2(R2). Amsterdam, pp. 1 – 34.
- Ivković, B., Marković, B., and Vladimirov, S. 2014. Development and Validation of RP-HPLC Method for Analysis of Multicomponent Cough-Cold Syrup Formulation. *Arh.farm.*, 64: 271 – 284.
- Jacobsen, D., and McMartin, K. E. (1986). Methanol and Ethylene Glycol Poisonings: Mechanism of Toxicity, Clinical Course, Diagnosis and Treatment. *Med. Toxicol.* 1: 309–334.
- Jones, D. 2008. *FASTtrack: Pharmaceuticals – Dosage Form and Design*. Pharmaceutical Press. London, UK.
- Leth, P.M., and Gregersen, M. 2005. Ethylene Glycol Poisoning. *Forensic Science International*. 155: 179 – 184. doi: 10.1016/j.forsciint.2004.11.012
- Li, H., Dong, J., Chen, W., Wang, S., Guo, H., Man, Y., et al. 2006. Measurement of Serum Total Glycerides and Free Glycerol by High-Performance Liquid Chromatography. *Journal of Lipid Research*. 47. (9): 2089 – 2096. doi: 10.1194/jlr.D600009-JLR200
- Liu, Y., Liu, Q., Wang, X., He, Z., Li, D., Guan, X., et al. 2018. Inhibition of Autophagy Attenuated Ethylene Glycol Induced Crystal Deposition and Renal Injury in a Rat Model of Nephrolithiasis. *Kidney and Blood Pressure Research*, 43(1): 246 – 255.
- Lopez-Sanchez, M., Dominguez-Vidal, A., Ayora-Canada, M.J., and Molina-Diaz, A. 2008. Assessment of Dentifrice Adulteration with Diethylene Glycol by Means of ATR-FTIR Spectroscopy and Chemometrics. *Analytical Chmica Acta*. 620: 113 – 119.
- Martinez, V.P.M., de la Garza, C.E.E., Hernandez-Garcia, A.G., and Perez, N.O. 2020. Detection and Quantification of Leached Ethylene Glycol in Biopharmaceuticals by RP-UHPLC. *Analytical and Bioanalytical Chemistry*. doi: 10.1007/s00216-020-02425-x
- Megarbane, B. 2010. Treatment of Patients with Ethylene Glycol or Methanol Poisoning: Focus on Fomepizole. *Open Access Emergency Medicine*. 2: 67 – 75. doi: 10.2147/OAEM.S5346
- Mikkelsen, S.R., and Corto´n, E. 2004. *Bioanalytical Chemistry. Chromatography of Biomolecules*. John Wiley & Sons, Inc., USA.
- Moosavi, S.M., and Ghassabian, S. 2018. Linearity of Calibration Curves for Analytical Methods: A Review of Criteria for Assessment of Method Reliability. *IntechOpen*. doi: 10.5772/intechopen.72932
- Nanco, C.R., Poklis, J.L., Hiler, M.M., Breland, A.B., Eissenberg, T., and Wolf, C.E. 2019. An Ultra-High-Pressure Liquid Chromatographic Tandem Mass Spectrometry Method for the Analysis of Benzoyl Ester Derivatized Glycols and Glycerol. *Journal of Analytical Toxicology*. 43: 720 – 725.


- Oehlke, J., Brudel, M., and Blasig, I.E. 1994. Benzoylation of Sugars, Polyols and Amino Acids in Biological Fluids for High Performance Liquid Chromatographic Analysis. *J. Chromatogr.* 655: 105 – 111.
- Rajpurohit, P., Boddu, G., and Malvi, A. 2022. Review Article: Diethylene Glycol Poisoning – Past to Present. *Journal of Pharmacovigilance and Drug Research.* 3. (4): 9 – 14. doi: 10.53411/jpadr.2022.3.4.3
- Rosano, T., Swift, T.A., Kranick, C.J., and Sikirica, M. 2009. Ethylene Glycol and Glycolic Acid in Postmortem Blood From Fatal Poisonings. *J Anal Toxicol.* 33. (8):508–13. doi : 10.1093/jat/33.8.508.
- Self, R. L. 2013. Direct Analysis in Real Time-Mass Spectrometry (DART-MS) for Rapid Qualitative Screening of Toxic Glycols in Glycerin-Containing Products. *Journal of Pharmaceutical and Biomedical Analysis.* 80: 155 – 158.
- Sharma, I., Khan, W., Parveen, R., Alam, M., Ahmad, I., Ansari, M.H.R., et al. 2017. Antiurolithiasis Activity of Bioactivity Guided Fraction of *Bergenia ligulate* Against Ethylene Glycol Induced Renal Calculi in Rat. *BioMed Research International*, 2017: 1969525
- Soe, J.W., Lee, J.H., Son, I.S., Kim, Y.J., Kim, D.Y., Hwang, Y., Chung, H.A., Choi, H.S., and Lim, S.D. 2012. Case Report: Acute Oxalate Nephropathy Caused by Ethylene Glycol Poisoning. *Kidney Research and Clinical Practice.* 31: 249 – 252. doi : 10.1016/j.krcp.2012.09.007
- Suhartati, T. 2017. *Dasar-Dasar Spektrofotometri UV-Vis dan Spektrometri Massa untuk Penentuan Struktur Senyawa Organik.* AURA. Bandar Lampung.
- USP Chapters <1086> Impurities in Drug Substance and Drug Products.
- Vollmer, P.A., Harty, D.C., Erickson, N.B., Balhon, A.C., and Dean, R.A. 1996. Serum Ethylene Glycol by High-Performance Liquid Chromatography. *J. Chromatogr.* 685: 370–374.
- Wang, Z. 2010. *Comprehensive Organic Name Reaction and Reagents.* John Wiley & Sons. USA.
- Wexler, P. 2024. *Encyclopedia of Toxicology.* Vol 1, 4th ed. Elsevier. US.
- Whalen, K., Field, C., and Radhakrishnan, R. 2019. *Pharmacology.* 7th ed. Wolters Kluwer. Philadelphia. p. 1668. Available as PDF file.
- World Health Organization (WHO). 2022. Medical Product Alert N°6/2022: Substandard (contaminated) paediatric medicines identified in WHO region of Africa. (Online). ([https://www.who.int/news/item/05-10-2022-medical-product-alert-n-6-2022-substandard-\(contaminated\)-paediatric-medicines](https://www.who.int/news/item/05-10-2022-medical-product-alert-n-6-2022-substandard-(contaminated)-paediatric-medicines), diakses 23 Februari 2023).
- Yao, H.H., and Porter, W.H.1996. Simultaneous Determination of Ethylene Glycol and Its Major Toxic Metabolite, Glycolic Acid, in Serum by Gas Chromatography. *Clin Chem.* 42. (2): 292–7.

Yulyarti, E., Rifai, Y., and Yulianty, R. 2018. Penetapan Kadar Parasetamol, Kafein dan Propifenazon secara Simultan dalam Sediaan Tablet dengan Metode KCKT. *Majalah Farmasi dan Farmakologi*. 22. (1): 1-4. <http://journal.unhas.ac.id/index.php/mff>


LAMPIRAN

Lampiran 1. Skema kerja umum

Pembuatan Kurva Baku



Analisis Residu EG dalam Sampel



Lampiran 2. Perhitungan

Perhitungan Linearitas

| Konsentrasi (ppm) | Waktu retensi (menit) | AUC |
|-------------------|-----------------------|---------|
| 1 | 2.805 | 223491 |
| 10 | 2.949 | 891804 |
| 20 | 2.998 | 1146924 |
| 30 | 3.069 | 1728335 |
| 40 | 3.086 | 2286412 |

Persamaan regresi:

$$y = a + bx$$

$$a = 233791$$

$$b = 50574$$

$$y = 233791 + 50574x$$

$$r = 0,992371$$

Perhitungan Batas Deteksi dan Batas Kuantitasi

| Konsentrasi (ppm) | AUC | y' | y-y' | y ² |
|-------------------|---------|------------------|--------|----------------|
| 1 | 223491 | 284365 | -60874 | 3705643876 |
| 10 | 891804 | 739531 | 152273 | 23187066529 |
| 20 | 1146924 | 1245271 | -98347 | 9672132409 |
| 30 | 1728335 | 1751011 | -22676 | 514200976 |
| 40 | 2286412 | 2256751 | 29661 | 879774921 |
| | | AVG | | 7591763742 |
| | | AVG(n) | | 1897940936 |
| | | Sy | | 43565,36394 |
| | | Batas deteksi | | 2,842680053 |
| | | Batas kuantitasi | | 8,61418198 |

$$\begin{aligned} \text{Batas deteksi} &= \frac{3,3 \text{ Sy}}{b} \\ &= \frac{3,3(43565,36394)}{50574} \\ &= 2,842680053 \text{ ppm} \end{aligned}$$

$$\text{Batas kuantitasi} = \frac{10 \text{ Sy}}{b}$$

$$= \frac{10(43565,36394)}{50574}$$

$$= 8,61418198 \text{ ppm}$$

Perhitungan Kadar EG pada Sampel

| Sampel | Waktu retensi (menit) | AUC | Kadar EG (ppm) | Kadar EG dalam Sirup (mg/mL) |
|--------|-----------------------|---------|------------------|------------------------------|
| A | 3,066 | 1733914 | 29,66194092 | 0,02966194092 |
| B | - | - | Tidak terdeteksi | Tidak terdeteksi |
| C | 2,832 | 343484 | 2,168960335 | 0,002168960335 |
| D | - | - | Tidak terdeteksi | Tidak terdeteksi |
| E | 2,941 | 183875 | -0,98698936 | -0,00098698936 |
| F | - | - | Tidak terdeteksi | Tidak terdeteksi |
| G | - | - | Tidak terdeteksi | Tidak terdeteksi |
| H | 2,907 | 613548 | 7,508937399 | 0,007508937399 |
| I | - | - | Tidak terdeteksi | Tidak terdeteksi |
| J | - | - | Tidak terdeteksi | Tidak terdeteksi |
| K | 2,911 | 509687 | 5,455293234 | 0,005455293234 |
| L | 3,027 | 1101220 | 17,15167873 | 0,01715167873 |
| M | - | - | Tidak terdeteksi | Tidak terdeteksi |
| N | 2,928 | 381006 | 2,910883062 | 0,002910883062 |
| O | - | - | Tidak terdeteksi | Tidak terdeteksi |

$$y = 233791 + 50574x$$

Sampel A

$$x = \frac{1733914 - 233791}{50574}$$

$$= 29,66194092 \text{ ppm}$$

$$= 29,66194092 \text{ mg/L}$$

$$= 0,02966194092 \text{ mg/mL}$$

Perhitungan Resolusi

$$R = \frac{2(t_{R2} - t_{R1})}{W_1 + W_2}$$

Resolusi untuk kromatogram 1 ppm

$$R = \frac{2(8,677 - 2,805)}{0,396 + 0,475}$$

$$R = 13,48335247$$

Perhitungan jumlah pelat teori (N)

$$N = 16 \left[\frac{t_R}{W} \right]^2$$

Nilai N untuk konsentrasi 1 ppm

$$N = 16 \left[\frac{2,805}{0,396} \right]^2$$

$$N = 802,7777778$$

Perhitungan HETP


$$H = \frac{L}{N}$$

Nilai HETP untuk konsentrasi 1 ppm


$$H = \frac{250 \text{ mm}}{802,7777778}$$

$$H = 0,311418685$$

Lampiran 3. Data Kromatogram




Gambar 5. Kromatogram baku konsentrasi 1 ppm



Gambar 6. Kromatogram baku konsentrasi 10 ppm

<Chromatogram>

mV




<Peak Table>

| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
|-------|-----------|---------|--------|--------|------|------|------|
| 1 | 2.998 | 1146924 | 96447 | 15.950 | | | |
| 2 | 8.675 | 1592017 | 163130 | 22.140 | | | |
| 3 | 17.005 | 4451689 | 243437 | 61.910 | | | |
| Total | | 7190630 | 503015 | | | | |

Gambar 7. Kromatogram baku konsentrasi 20 ppm

<Chromatogram>

mV




<Peak Table>

| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
|-------|-----------|----------|--------|--------|------|------|------|
| 1 | 3.069 | 1728335 | 109413 | 12.936 | | | |
| 2 | 8.637 | 2329081 | 238899 | 17.433 | | | |
| 3 | 16.692 | 9302889 | 516368 | 69.631 | | | |
| Total | | 13360304 | 864680 | | | | |

Gambar 8. Kromatogram baku konsentrasi 30 ppm

<Chromatogram>

mV




<Peak Table>

| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
|-------|-----------|---------|--------|--------|------|------|------|
| 1 | 3.086 | 2286412 | 202514 | 28.773 | | | |
| 2 | 8.680 | 1409781 | 143797 | 17.741 | | | |
| 3 | 17.035 | 4250058 | 231608 | 53.485 | | | |
| Total | | 7946252 | 577919 | | | | |

Gambar 9. Kromatogram baku konsentrasi 40 ppm

<Chromatogram>

mV




<Peak Table>

| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
|-------|-----------|---------|--------|--------|------|------|------|
| 1 | 3.066 | 1733914 | 122699 | 20.605 | | SV | |
| 2 | 6.340 | 178560 | 18179 | 2.122 | | V | |
| 3 | 8.365 | 398235 | 39239 | 4.732 | | | |
| 4 | 8.731 | 1291889 | 132933 | 15.352 | | V | |
| 5 | 9.224 | 135319 | 9577 | 1.608 | | V | |
| 6 | 11.525 | 112933 | 6508 | 1.342 | | | |
| 7 | 12.352 | 275626 | 13599 | 3.275 | | SV | |
| 8 | 17.614 | 3648372 | 163023 | 43.355 | | | |
| 9 | 19.231 | 640206 | 27014 | 7.608 | | V | |
| Total | | 8415053 | 532774 | | | | |

Gambar 10. Kromatogram sampel A

<Chromatogram>

mV




<Peak Table>

| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
|-------|-----------|----------|---------|--------|------|------|------|
| 1 | 3.158 | 5365449 | 368925 | 25.830 | | | |
| 2 | 7.680 | 1797875 | 58645 | 8.655 | | | |
| 3 | 8.726 | 7813457 | 802109 | 37.615 | | S | |
| 4 | 12.301 | 3479671 | 94047 | 16.752 | | | |
| 5 | 17.187 | 2315591 | 118163 | 11.148 | | V | |
| Total | | 20772044 | 1441888 | | | | |

Gambar 11. Kromatogram sampel B

<Chromatogram>

mV




<Peak Table>

| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
|-------|-----------|----------|---------|--------|------|------|------|
| 1 | 2.832 | 343484 | 32702 | 1.255 | | V | |
| 2 | 3.657 | 1612297 | 55349 | 5.890 | | V | |
| 3 | 4.506 | 445620 | 16448 | 1.628 | | V | |
| 4 | 6.332 | 749465 | 76312 | 2.738 | | V | |
| 5 | 7.649 | 1097154 | 35413 | 4.008 | | V | |
| 6 | 11.622 | 21056771 | 1623440 | 76.926 | | V | |
| 7 | 12.322 | 409671 | 20962 | 1.497 | | V | |
| 8 | 17.550 | 538857 | 25844 | 1.969 | | | |
| 9 | 19.340 | 1119561 | 21212 | 4.090 | | | |
| Total | | 27372880 | 1907682 | | | | |

Gambar 12. Kromatogram sampel C

<Chromatogram>




<Peak Table>

| Detector A 232nm | | | | | | | |
|------------------|-----------|----------|--------|--------|------|------|------|
| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
| 1 | 2.689 | 40447 | 5591 | 0.181 | | V | |
| 2 | 4.442 | 43919 | 2389 | 0.196 | | V | |
| 3 | 5.917 | 2838084 | 82733 | 12.684 | | V | |
| 4 | 6.123 | 1823998 | 80615 | 8.152 | | V | |
| 5 | 7.591 | 12867721 | 315050 | 57.510 | | SV | |
| 6 | 8.722 | 827085 | 83403 | 3.697 | | V | |
| 7 | 9.199 | 196571 | 12084 | 0.879 | | V | |
| 8 | 10.485 | 75250 | 5611 | 0.336 | | | |
| 9 | 11.608 | 1016919 | 57357 | 4.545 | | | |
| 10 | 12.328 | 75597 | 4131 | 0.338 | | V | |
| 11 | 13.031 | 42372 | 717 | 0.189 | | V | |
| 12 | 17.563 | 470222 | 23259 | 2.102 | | | |
| 13 | 19.261 | 2056401 | 38055 | 9.191 | | | |
| Total | | 22374586 | 710997 | | | | |

Gambar 13. Kromatogram sampel D

<Chromatogram>




<Peak Table>

| Detector A 232nm | | | | | | | |
|------------------|-----------|----------|---------|--------|------|------|------|
| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
| 1 | 2.941 | 183875 | 25736 | 0.318 | | | |
| 2 | 3.659 | 8625281 | 165419 | 14.938 | | V | |
| 3 | 4.621 | 10084255 | 146940 | 17.465 | | V | |
| 4 | 6.125 | 314349 | 28047 | 0.544 | | V | |
| 5 | 6.325 | 142547 | 14263 | 0.247 | | V | |
| 6 | 8.716 | 6333196 | 653387 | 10.969 | | V | |
| 7 | 11.614 | 27765016 | 2204590 | 48.087 | | S | |
| 8 | 17.146 | 4290124 | 89361 | 7.430 | | V | |
| Total | | 57738642 | 3327743 | | | | |

Gambar 14. Kromatogram sampel E

<Chromatogram>

mV



<Peak Table>


Detector A 232nm

| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
|-------|-----------|----------|---------|--------|------|------|------|
| 1 | 2.783 | 208219 | 22144 | 0.486 | | V | |
| 2 | 3.599 | 6458334 | 205133 | 15.071 | | V | |
| 3 | 4.493 | 1483101 | 25416 | 3.461 | | V | |
| 4 | 6.333 | 1268903 | 122933 | 2.961 | | SV | |
| 5 | 7.601 | 1061975 | 36605 | 2.478 | | V | |
| 6 | 8.086 | 560595 | 48828 | 1.308 | | V | |
| 7 | 8.342 | 94961 | 9481 | 0.222 | | V | |
| 8 | 8.718 | 482869 | 48523 | 1.127 | | V | |
| 9 | 11.617 | 30270628 | 2386653 | 70.637 | | SV | |
| 10 | 16.567 | 672273 | 30301 | 1.569 | | V | |
| 11 | 17.529 | 292034 | 13932 | 0.681 | | V | |
| Total | | 42853893 | 2949949 | | | | |

Gambar 15. Kromatogram sampel F

<Chromatogram>

mV



<Peak Table>


Detector A 232nm

| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
|-------|-----------|----------|---------|--------|------|------|------|
| 1 | 1.221 | 123082 | 4657 | 0.210 | | | |
| 2 | 3.161 | 4971613 | 391564 | 8.483 | | | |
| 3 | 3.558 | 3787446 | 105983 | 6.462 | | V | |
| 4 | 4.264 | 2106912 | 68389 | 3.595 | | SV | |
| 5 | 6.141 | 5979328 | 688370 | 10.202 | | V | |
| 6 | 7.564 | 21303931 | 700376 | 36.349 | | SV | |
| 7 | 8.359 | 53729 | 7239 | 0.092 | | T | |
| 8 | 8.713 | 5784914 | 595471 | 9.870 | | V | |
| 9 | 9.073 | 789832 | 65598 | 1.348 | | SV | |
| 10 | 10.882 | 25368 | 793 | 0.043 | | | |
| 11 | 11.621 | 55781 | 4039 | 0.095 | | V | |
| 12 | 12.082 | 93073 | 2309 | 0.159 | | V | |
| 13 | 13.740 | 56563 | 3268 | 0.097 | | V | |
| 14 | 14.945 | 2234819 | 123167 | 3.813 | | V | |
| 15 | 17.144 | 1963739 | 108669 | 3.351 | | | |
| 16 | 17.520 | 2064255 | 104477 | 3.522 | | V | |
| 17 | 19.103 | 7215576 | 123031 | 12.311 | | | |
| Total | | 58609963 | 3097401 | | | | |

Gambar 16. Kromatogram sampel G

<Chromatogram>

mV




<Peak Table>

| Detector A 232nm | | | | | | | |
|------------------|-----------|---------|--------|--------|------|------|------|
| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
| 1 | 2.907 | 613548 | 50700 | 17.150 | | V | |
| 2 | 4.260 | 66244 | 10122 | 1.852 | | V | |
| 3 | 6.135 | 266258 | 35737 | 7.442 | | V | |
| 4 | 6.319 | 633305 | 65618 | 17.702 | | V | |
| 5 | 7.493 | 77655 | 2581 | 2.171 | | V | |
| 6 | 8.059 | 1162320 | 105509 | 32.489 | | V | |
| 7 | 8.692 | 340589 | 34680 | 9.520 | | V | |
| 8 | 9.084 | 87976 | 5304 | 2.459 | | V | |
| 9 | 14.208 | 66749 | 1759 | 1.866 | | | |
| 10 | 14.879 | 160023 | 5046 | 4.473 | | V | |
| 11 | 17.430 | 102954 | 4933 | 2.878 | | | |
| Total | | 3577621 | 321988 | | | | |

Gambar 17. Kromatogram sampel H

<Chromatogram>

mV




<Peak Table>

| Detector A 232nm | | | | | | | |
|------------------|-----------|---------|--------|--------|------|------|------|
| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
| 1 | 2.662 | 32341 | 3967 | 0.369 | | V | |
| 2 | 3.423 | 833695 | 29297 | 9.502 | | V | |
| 3 | 4.241 | 25984 | 1177 | 0.296 | | V | |
| 4 | 4.603 | 11541 | 1385 | 0.132 | | V | |
| 5 | 5.236 | 25134 | 974 | 0.286 | | V | |
| 6 | 6.133 | 71368 | 9152 | 0.813 | | V | |
| 7 | 6.316 | 110865 | 11768 | 1.264 | | V | |
| 8 | 7.396 | 6627371 | 219069 | 75.534 | | SV | |
| 9 | 8.689 | 111684 | 11569 | 1.273 | | T | |
| 10 | 10.422 | 107071 | 8193 | 1.220 | | S | |
| 11 | 11.423 | 767367 | 49724 | 8.746 | | SV | |
| 12 | 18.719 | 49551 | 796 | 0.565 | | | |
| Total | | 8773971 | 347070 | | | | |

Gambar 18. Kromatogram sampel I

<Chromatogram>

mV




<Peak Table>

| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
|-------|-----------|-----------|---------|--------|------|------|------|
| 1 | 3.149 | 7911746 | 493198 | 6.909 | | | |
| 2 | 4.008 | 699270 | 56299 | 0.611 | | V | |
| 3 | 4.260 | 992236 | 60068 | 0.867 | | V | |
| 4 | 4.602 | 1048113 | 66818 | 0.915 | | V | |
| 5 | 5.779 | 37226518 | 3999587 | 32.510 | | V | |
| 6 | 6.312 | 2189262 | 209528 | 1.912 | | V | |
| 7 | 6.698 | 52709 | 5123 | 0.046 | | V | |
| 8 | 7.375 | 3636023 | 117928 | 3.175 | | SV | |
| 9 | 8.685 | 19856008 | 2025504 | 17.340 | | V | |
| 10 | 9.173 | 2489736 | 209441 | 2.174 | | V | |
| 11 | 11.566 | 414173 | 27981 | 0.362 | | V | |
| 12 | 12.325 | 304284 | 5560 | 0.266 | | V | |
| 13 | 14.863 | 185115 | 9908 | 0.162 | | | |
| 14 | 16.313 | 93061 | 5349 | 0.081 | | V | |
| 15 | 17.442 | 37409285 | 1532650 | 32.670 | | SV | |
| Total | | 114507539 | 8824941 | | | | |

Gambar 19. Kromatogram sampel J

<Chromatogram>




<Peak Table>

| Detector A 232nm | | Area | Height | Conc. | Unit | Mark | Name |
|------------------|-----------|----------|---------|--------|------|------|------|
| Peak# | Ret. Time | | | | | | |
| 1 | 2.604 | 27853 | 2145 | 0.206 | | | |
| 2 | 2.911 | 509687 | 46287 | 3.771 | | V | |
| 3 | 4.260 | 13008 | 1928 | 0.096 | | V | |
| 4 | 5.415 | 11752 | 1432 | 0.087 | | V | |
| 5 | 6.308 | 4540209 | 482924 | 33.595 | | S | |
| 6 | 7.415 | 90748 | 9553 | 0.671 | | T | |
| 7 | 8.042 | 167009 | 15852 | 1.236 | | T | |
| 8 | 8.669 | 704406 | 72758 | 5.212 | | V | |
| 9 | 9.143 | 284540 | 24429 | 2.105 | | V | |
| 10 | 11.536 | 4341196 | 347470 | 32.122 | | | |
| 11 | 12.196 | 80257 | 5182 | 0.594 | | V | |
| 12 | 12.804 | 12467 | 818 | 0.092 | | V | |
| 13 | 17.338 | 2731531 | 127118 | 20.212 | | | |
| Total | | 13514662 | 1137896 | | | | |

Gambar 20. Kromatogram sampel K

<Chromatogram>




<Peak Table>

| Detector A 232nm | | Area | Height | Conc. | Unit | Mark | Name |
|------------------|-----------|----------|---------|--------|------|------|------|
| Peak# | Ret. Time | | | | | | |
| 1 | 1.111 | 28065 | 984 | 0.091 | | | |
| 2 | 3.027 | 1101220 | 79130 | 3.561 | | | |
| 3 | 3.256 | 318703 | 17400 | 1.030 | | V | |
| 4 | 4.106 | 98799 | 3447 | 0.319 | | SV | |
| 5 | 6.065 | 51743 | 6466 | 0.167 | | V | |
| 6 | 6.299 | 540447 | 56830 | 1.747 | | V | |
| 7 | 7.169 | 25147153 | 838052 | 81.310 | | SV | |
| 8 | 8.025 | 152619 | 16285 | 0.493 | | T | |
| 9 | 8.657 | 1706749 | 175220 | 5.519 | | SV | |
| 10 | 11.532 | 531923 | 16381 | 1.720 | | | |
| 11 | 12.153 | 69328 | 4493 | 0.224 | | V | |
| 12 | 14.682 | 19156 | 332 | 0.062 | | | |
| 13 | 16.924 | 1161516 | 64227 | 3.756 | | | |
| Total | | 30927420 | 1279247 | | | | |

Gambar 21. Kromatogram sampel L

<Chromatogram>

mV



<Peak Table>


Detector A 232nm

| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
|-------|-----------|----------|---------|--------|------|------|------|
| 1 | 3.100 | 2627545 | 186784 | 15.008 | | | |
| 2 | 6.294 | 2990635 | 316511 | 17.082 | | V | |
| 3 | 7.097 | 1049905 | 33130 | 5.997 | | V | |
| 4 | 8.020 | 328358 | 27898 | 1.876 | | V | |
| 5 | 8.652 | 2649319 | 206161 | 15.132 | | SV | |
| 6 | 11.593 | 112561 | 5916 | 0.643 | | T | |
| 7 | 12.154 | 193747 | 13146 | 1.107 | | TV | |
| 8 | 13.942 | 131595 | 2194 | 0.752 | | TV | |
| 9 | 16.026 | 292012 | 7679 | 1.668 | | TV | |
| 10 | 16.904 | 5780650 | 308431 | 33.018 | | V | |
| 11 | 17.782 | 1351194 | 24155 | 7.718 | | V | |
| Total | | 17507521 | 1132007 | | | | |

Gambar 22. Kromatogram sampel M

<Chromatogram>

mV



<Peak Table>


Detector A 232nm

| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
|-------|-----------|----------|---------|--------|------|------|------|
| 1 | 1.085 | 83832 | 3183 | 0.387 | | | |
| 2 | 2.928 | 381006 | 55113 | 1.757 | | | |
| 3 | 3.319 | 2125835 | 65693 | 9.805 | | V | |
| 4 | 3.876 | 1655808 | 42817 | 7.637 | | V | |
| 5 | 5.239 | 1660589 | 23159 | 7.659 | | V | |
| 6 | 6.293 | 925820 | 92587 | 4.270 | | V | |
| 7 | 7.071 | 3627051 | 116441 | 16.729 | | SV | |
| 8 | 8.647 | 3989971 | 412625 | 18.403 | | | |
| 9 | 9.111 | 120632 | 9778 | 0.556 | | V | |
| 10 | 11.494 | 27650 | 1790 | 0.128 | | V | |
| 11 | 16.888 | 7000060 | 380250 | 32.286 | | | |
| 12 | 18.827 | 83067 | 3653 | 0.383 | | | |
| Total | | 21681323 | 1207090 | | | | |

Gambar 23. Kromatogram sampel N

<Chromatogram>

mV



<Peak Table>

Detector A 232nm

| Peak# | Ret. Time | Area | Height | Conc. | Unit | Mark | Name |
|-------|-----------|----------|---------|--------|------|------|------|
| 1 | 2.719 | 83320 | 10530 | 0.320 | | V | |
| 2 | 5.015 | 697006 | 10623 | 2.674 | | | |
| 3 | 5.397 | 288391 | 11734 | 1.106 | | V | |
| 4 | 6.291 | 14463097 | 1533801 | 55.485 | | V | |
| 5 | 6.628 | 669653 | 66261 | 2.569 | | V | |
| 6 | 7.058 | 998272 | 34250 | 3.830 | | SV | |
| 7 | 8.224 | 102521 | 8434 | 0.393 | | V | |
| 8 | 8.641 | 414210 | 35269 | 1.589 | | V | |
| 9 | 9.102 | 305937 | 24344 | 1.174 | | V | |
| 10 | 11.322 | 203775 | 12031 | 0.782 | | V | |
| 11 | 13.857 | 212523 | 4279 | 0.815 | | V | |
| 12 | 14.617 | 83804 | 2993 | 0.321 | | V | |
| 13 | 15.158 | 112974 | 2527 | 0.433 | | V | |
| 14 | 17.224 | 3303644 | 114312 | 12.674 | | V | |
| 15 | 18.804 | 4127779 | 177296 | 15.835 | | V | |
| Total | | 26066906 | 2048685 | | | | |

Gambar 24. Kromatogram sampel O

Lampiran 4. Dokumentasi Penelitian**Gambar 25. Menimbang bahan****Gambar 26. Penyiapan larutan baku****Gambar 27. Sampel sirup flu dan batuk anak**



Gambar 28. Penyaringan larutan

CURRICULUM VITAE

A. Data pribadi

1. Nama : Aliyah Sukma
2. Tempat, tgl. Lahir : Siwa, 22 Mei 2001
3. Alamat : Jl. Politeknik, Tamalanrea
4. Kewarganegaraan : Warga Negara Indonesia

B. Riwayat Pendidikan

1. Tamat SMP tahun 2016 di MTsN 1 Wajo
2. Tamat SMA tahun 2019 di SMAN 6 Wajo

C. Pekerjaan dan Riwayat Pekerjaan

- Jenis pekerjaan : -
- NIP atau identitas lain (NIK) : -
- Pangkat/Jabatan : -

D. Karya ilmiah yang telah dipublikasikan (misalnya pada jurnal):

-

E. Makalah pada Seminar/Konferensi Ilmiah Nasional dan Internasional

-